

EasyScript® One-Step RT-PCR SuperMix

Cat. No. AE411

Storage: at -20°C for two years

Description

One-Step RT-PCR combines the first-strand cDNA synthesis with PCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene-specific primers can be used for One-Step RT-PCR. *EasyScript*® RT and *TransTaq*® HiFi DNA Polymerase are used in the kit.

Highlights

Amplification of fragment up to 4 kb.

Applications

Multiple copy gene detection.

Kit Contents

Component	AE411-02 (200 rxns)
<i>EasyScript</i> ® One-Step Enzyme Mix	80 µl
2×ES One-Step Reaction Mix	2×1 ml
RNase-free Water	2×1 ml

Reaction Components

Component	Volume	Final Concentration
RNA Template	1 pg~1 µg	as required
Forward GSP(10 µM)	0.4 µl	0.2 µM
Reverse GSP(10 µM)	0.4 µl	0.2 µM
2×ES One-Step Reaction Mix	10 µl	1×
<i>EasyScript</i> ® One-Step Enzyme Mix	0.4 µl	-
RNase-free Water	Variable	-
Total volume	20 µl	-

Thermal cycling conditions

45°C	15-30 min	} 35-40 cycles
94°C	2-5 min	
94°C	30 sec	
50-60°C	30 sec	
72°C	1-2 kb/min	
72°C	5-10 min	

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

