

RNase-free Water

Please read the manual carefully before use.

Cat. No. GI201

Storage: at 15°C-30°C for two years

Description

RNase-free Water is prepared from deionized water incubated with 0.01% DEPC, then autoclaved to remove residual DEPC. It is suitable for various RNA-related experiments.

Notes

This product has been autoclaved and no longer contains DEPC. Please avoid RNase contamination during use.

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email custserv@transgenbiotech.com

