

## Penicillin-Streptomycin (100×)

**Cat. No.** FG101

**Storage:** at -18°C or below for one year

**Description**

Penicillin-Streptomycin (100×) contains 10 kU/ml of penicillin and 10 mg/ml of streptomycin. The solution has been filter-sterilized. It can be used for cell culture at a final concentration of 1x.

**Note**

Aliquots after thawing; avoiding repeated freezing and thawing; store at 2-8°C for two weeks.

**For research use only, not for clinical diagnosis**

Service telephone +86-10-57815020

Service email [custserv@transgenbiotech.com](mailto:custserv@transgenbiotech.com)

